

# C-Chip

## INSTRUCTIONS

Disposable hemocytometer

---

System Neubauer Improved

---

DHC-N01

---



Website: <http://www.digital-bio.com>

E-mail: [sales@digital-bio.com](mailto:sales@digital-bio.com)

NanoEnTek Inc.

Worldwide:

12F, Ace High-end Tower, 235-2 Guro-3dong, Guro-gu, Seoul, 152-711, Korea

Tel: +82-2-6220-7911 Fax: +82-2-6220-7721

USA:

63 South St., Suite 250, Hopkinton, MA 01748, USA

Tel: +1-508-544-9809 Fax: +1-508-544-1554



731-30, Jiwol-ri, Paltan-Myeon, Hwasung, Gyeonggi, 445-911, Korea

---

All the materials in this manual are protected by Korean and international copyright laws.

This manual may not be reproduced, translated, published or distributed by any means without the permission of the copyright owner.

The information in this manual is described as correctly as possible, but may be changed without prior consent or notification.

Documentation : DHC-N01(Rev NI.00.00-02/2008)

Copyright © 2008, by NanoEnTek Inc.

All rights reserved. Published in Korea.

## Unpacking

When you receive the C-Chip (DHC-N01) for the first time, you will find the following components in your package.

Disposable hemocytometer  
Instruction manual

## Safety Precautions

For analyzing hazardous or potential infectious materials:

Take necessary precautions

Handle with care

Dispose in an appropriate way

Long exposure to solvents will cause the slide to warp. Xylene and toluene based mounting media should be avoided. Glycerol, gelatin, and other aqueous-based media are recommended.

## Safety Symbols:

The safety symbols on the C-Chip (DHC-N01) are intended to inform you. For potential danger or a particular caution. Before use, please read and consult the guide for the symbols and their meanings.



Batch Code (Lot Number)



Expire Date



Do not Re-use



A Manufacturer



See Instruction for Use


**NOTE:** The C-Chip (DHC-N01) is for **single use** only. **Do not reuse.** It should be used immediately after unsealing. The warranty on the C-Chip included in the conditions of supply is valid for 24 months from the date of manufacturing. The **expiration date** is indicated on the front of outer box.

## Introduction


The C-Chip (DHC-N01) is a disposable plastic hemocytometer used for manual cell counting. It consists of surface-patterned two enclosed chambers with two ports for sample injection (Fig. 2).

The DHC-N01 grid pattern is exactly same as the Neubauer Improved. It consists of 9 large squares, each measuring 1 x 1 mm, and the depth of the chamber is 0.1 mm. Each square has a total volume of  $0.1 \text{ mm}^3$  or  $10^{-4} \text{ cm}^3$  (Fig. 1).

The central square is divided into 25 small squares with triple lines and four corner squares are divided into 16 small squares.



*Figure 1. Grid pattern*



*Figure 2. DHC-N01*

## Counting with C-Chip

### A. General Methods

1. Mix the samples well.
2. Load 10  $\mu\text{l}$  of sample into the sample injection area in Fig. 2, so that it fills the chamber by capillary action.
3. Count the cells under the microscope.

$$\text{Cells per ml} = \text{Average count per square} \times \text{dilution factor} \times \text{volume factor}$$

### B. Mammalian cell counting

1. Treat the cell samples with trypsin-EDTA.
2. Carefully remove the supernatant with a pipette tip without disturbing the pellet.
3. Add an appropriate volume of growth media or PBS to dilute to a final concentration of  $5 \times 10^3$  cells/ml to  $5 \times 10^6$  cells per ml.
4. Thoroughly resuspend the cell pellet with a pipette.
5. Check visually if there are any cell clumps or agglomerates.
6. Load 10  $\mu\text{l}$  of sample into the sample injection area in Fig. 2.
7. Count the cells in 5 large squares.

Cells per ml =

$$\frac{\text{cells in 5 large squares}}{5} \times \text{dilution factor} \times 10^4 (\text{volume factor})$$

### C. Erythrocyte counting (1:200 dilution)

1. Dilute blood using accepted laboratory methods.
2. Load 10 µl of diluted sample into the sample injection area in Fig. 2.
3. Count the erythrocytes in the 5 small squares (four small corner squares and one small middle square) of the large center square.

RBCs per ml =

Cells in 5 small squares  $\times$  5  $\times$  200 (dilution factor)  $\times$   $10^4$  (volume factor)

### D. Leukocyte counting (1:20 dilution)

1. Dilute blood using accepted laboratory methods.
2. Load 10 µl of diluted sample into the sample injection area in Fig. 2.
3. Count the leukocytes in the 4 large corner squares.

WBCs per ml =

$$\frac{\text{cells in 4 corner squares}}{4} \times 20(\text{dilution factor}) \times 10^4 (\text{volume factor})$$

### Trouble shooting

In case of poor visibility results.

- ✓ Carefully load samples into the C-Chip in order to prevent the introduction of air bubbles.
- ✓ Observe after removing the dust from samples.
- ✓ Adjust the focus of the microscope.