

# Performance Evaluation of CAR-T Cell Source Isolated and Activated T Cells Using the ADAM-CellIT

Kyuheon Lee, Jonghyeon Oh, SangHwa Ahn, Hyung Seop Lee, Jaeseung Jeng.

NanoEnTek, Seoul, Republic of Korea

## Introduction


- Chimeric antigen receptor (CAR) T cell one of the immune cell therapy products, which is used for disease therapy such as cancer by overexpressing specific gene on T cells.
- Recently, CAR-T cells are being touted as the next generation cell to increase treatment efficiency<sup>1</sup>.
- For that reason, the importance of monitoring of the intermediate products has come into the limelight<sup>2</sup>.
- However, GMP facilities lack proper equipment which validates intermediate products, isolated or activated T cell, during CAR-T cell manufacturing processes<sup>3</sup>.
- ADAM-CellIT has been developed with the purpose of monitoring the quantity and viability of intermediate products in CAR-T cell manufacturing processes and also it is an optimized product for accurately and quickly counting isolated and activated T cells.

## Methods and Procedures


- The precision, linearity, and method comparison of ADAM-CellIT were evaluated in accordance with CLSI guidelines EP05-A3, EP06-A, and EP09-A3.
- The isolated and activated T cells were made from human blood and were measured both with ADAM-CellIT and by a comparative assay.
- The comparative assay was performed according to BD Leucocount™ Kit (BD biosciences; 340523) protocol using flow cytometry.

## Optimization Workflow for Producing CAR-T

### Principle of Viability Measurement (PI-staining Method)



### QC Platform for Producing CAR-T



- It is easy to monitor all different steps of the purification, expansion, and formulation of CAR-T cells using the ADAM-MC2 and ADAM-CellIT to ensure precise and reliable results. ADAM-MC2 and ADAM-CellIT can be used for cGMP, process control and quality control of CAR-T cell manufacturing.


## Conclusion

- Data indicates that the newly developed ADAM-CellIT assay exhibits reliable performance. Consequently, we expect that ADAM-CellIT will be a useful equipment to manage the monitoring of the quality of intermediate products during CAR-T cell manufacturing processes.

## Results

- The Coefficients of variation of precision within each measuring series did not exceed 10% (A, D).
- The analytical measurement range of the assays were  $5 \times 10^4 \sim 4 \times 10^6$  cell/mL with ordinary least squares regression fit of  $y = 1.0196x - 27468$  ( $r^2=0.9962$ ) (B).
- In the method comparison studies with BD Leucocount™ Kit, the correlation coefficient ( $r$ ) was 0.994, and the slopes/intercepts were 1.008 (95% CI= 0.9836 to 1.029)/-114.2(95% CI= -45499 to 26091) by Passing-Bablok regression fit (C).

### Performance Evaluation in Isolated T Cells



The profiles of T cells which was used for performance evaluation - Flow cytometric analysis of CD3 expression on unsorted (WBC; left panel) or sorted (Isolated T cell; right panel) human peripheral blood lymphocytes

#### A Precision

CLSI guideline EP05-A3

Level = Low

Level	Replicates	Mean	SD	CV
Low	20	486,336	24,181	4.97%

Level = Medium


Level	Replicates	Mean	SD	CV
Medium	20	1,446,250	44,073	3.05%

Level = High

Level	Replicates	Mean	SD	CV
High	20	2,905,800	48,043	1.65%


The sample that low, medium, and high concentration of isolated T cells were counted with ADAM-CellIT

#### B Linearity



Comparison between flow cytometry and ADAM-CellIT in Isolated T cells.


#### C



Correlation of T cell counting counting between flow cytometry and ADAM-CellIT in Isolated T cells.

Sample size	79
$y = -114.2 + 1.008x$	
<b>Systematic differences</b>	
Intercept A	-114.2
95% CI	-45499 to 26091
<b>Proportional differences</b>	
Slope B	1.008
95% CI	0.9836 to 1.029

### Performance Evaluation in Activated T Cells



The phenotypes of activated T cells which was used for performance evaluation - Flow cytometric analysis of CD25 expression on TCR stimulated (Activated T cell; bottom panel) or unstimulated (control T cell; upper panel) T cells.

#### D Precision

CLSI guideline EP05-A3

Level = Low

Level	Replicates	Mean	SD	CV
Low	20	514,470	48,655	9.46%

Level = Medium


Level	Replicates	Mean	SD	CV
Medium	20	1,436,925	92,938	6.47%

Level = High

Level	Replicates	Mean	SD	CV
High	20	2,887,715	145,705	5.05%

The sample that low, medium, and high concentration of activated T cells were counted with ADAM-CellIT

#### E



Comparison of viability between flow cytometry and ADAM-CellIT in Activated T cells

## References

1. Riddell SR. *Cancer J.* 20, 141-144 (2014).
2. Yonghong Li, *Engineering.* 5, 122-131 (2019).
3. Gee AP. *Best Pract Res Clin Haematol.* 31, 126-134 (2018).