


# Disposable Sampler (AMIES)



## Features

- Simplify specimen collection, transport and processing.
- Fewer samples need to be collected offering a more comfortable experience.
- Better preservation ability: the colony number of samples was reduced by no more than  $1.5\log_{10}$  after preservation for 48h.

## Application



Automation



Rapid antigen tests



Molecular testing



Gram stains


NEST disposable sampler (AMIES) is a collection and transport system which is used for aerobic, anaerobic and fastidious bacteria. In addition, NEST disposable sampler (AMIES) can be used for nucleic acid test of bacteria.

## Storage and Transport Temperature Requirement

**Before using:** AMIES solution can be stored in room temperature (5 to 25 degree centigrade) for twelve months.

**After using:** sample could maintain viability for up to 48 hours at room or refrigerator temperature. (Neisseria gonorrhoeae survival at 24 hours per CLSI standard).

## Quick Guide



### Collect

Collect the sample with the swab. Avoid touching the swab applicator below the breakpoint.



#### Oropharyngeal Swab & Capture Cap Swab:

- Gently insert swab into back of the throat and tonsillar area.
- Rub the swab over both tonsillar pillars and the posterior oropharynx, avoid touching the tongue, teeth and gums.



#### Nasopharyngeal Swab:

- With the patient seated, if possible, tilt the head back 70 degrees.
- Gently insert swab into the nostril along the septum floor of the nose extending straight back until the posterior nasopharynx is reached (distance from nostrils to external opening of ear).
- Rotate the swab several times while the swab is in contact with the nasopharyngeal wall.



### Snap

- Remove cap from the vial and insert the swab to the bottom of the vial.
- Holding the vial away from face, grasp the end of swab shaft and bend at a 180 degree angle to break at the breakpoint.
- Screw the cap on tightly to prevent leakage.




### Mark

Write information which you want on the vial body.


## Capture Cap



| Cat.No. | Description   | /Pk | /Cs |
|---------|---|-----|-----|
| 202451  | Disposable Sampler, 5 mL Vial with 1.2 mL ATM, with Individually Wrapped and Sterile Capture Cap Swabs, 10 Vials/pk, 50 Vials + 50 Swabs/cs | 10  | 50  |
| 202452  | Disposable Sampler, 5 mL Vial with 1.2 mL ATM, with Individually Wrapped and Sterile Capture Cap Swabs, 1 Vial + 1 Swab/pk, 100 pk/cs       | 1   | 100 |
| 202491  | Disposable Sampler, 5 mL Vial (empty vial), with Individually Wrapped and Sterile Capture Cap Swabs, 1 Vial + 1 Swab/pk, 100 pk/cs          | 1   | 100 |



5 mL Capture Cap




Capture Cap Swabs


## Normal Cap




| Cat.No. | Description   | /Pk | /Cs |
|---------|---|-----|-----|
| 202051  | Disposable Sampler, 5 mL Vial with 1.2 mL ATM, with Individually Wrapped and Sterile Oropharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs  | 10  | 50  |
| 202055  | Disposable Sampler, 5 mL Vial with 1.2 mL ATM, with Individually Wrapped and Sterile Nasopharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs | 10  | 50  |
| 202081  | Disposable Sampler, 10 mL Vial with 3 mL ATM, with Individually Wrapped and Sterile Oropharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs   | 10  | 50  |
| 202085  | Disposable Sampler, 10 mL Vial with 3 mL ATM, with Individually Wrapped and Sterile Nasopharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs  | 10  | 50  |
| 202052  | Disposable Sampler, 5 mL Vial with 1.2 mL ATM, with Individually Wrapped and Sterile Oropharyngeal Swabs, 1 Vial + 1 Oral Swab/pk, 100 pk/cs  | 1   | 100 |
| 202054  | Disposable Sampler, 5 mL Vial with 1.2 mL ATM, with Individually Wrapped and Sterile Nasopharyngeal Swabs, 1 Vial + 1 NP Swab/pk, 100 pk/cs   | 1   | 100 |
| 202056  | Disposable Sampler, 10 mL Vial with 3 mL ATM, with Individually Wrapped and Sterile Oropharyngeal Swabs, 1 Vial + 1 Oral Swab/pk, 100 pk/cs   | 1   | 100 |
| 202057  | Disposable Sampler, 10 mL Vial with 3 mL ATM, with Individually Wrapped and Sterile Nasopharyngeal Swabs, 1 Vial + 1 NP Swab/pk, 100 pk/cs    | 1   | 100 |



10 mL Normal Cap



5 mL Normal Cap



Oropharyngeal Swabs