

## ORBITAL-GENIE

## **Benchtop Orbital Shaker**



- O Choice of THREE different large footprint units depending on preferred platform (18 x 18in).
- Variable shaking speed with orbit REVERSING.
- Fully programmable time, reversing time, acceleration/deceleration rates, and variable RUN/STOP time.
- Create any mixing action from a smooth vortex to sloshing aeration (or anything in between).
- Programmable ramping feature ensures a smooth ramp-up to desired speed avoiding potential spills.
- Large independent LED's for speed and time allowing operator to view all settings at once.
- May be used in cold room or incubator and CO2 environments.



- 1. Orbital-Genie With Universal Ratcheting Clamps provides the flexibility for using a mix of different size and shape vessels on a single platform. Unit is supplied with two ratcheting and one non-ratcheting cushioned clamps. Extra ratcheting clamps are available.
- 2. Orbital-Genie With Adhering Mat provides a flat platform with adhering rubber mat perfect for use with open containers. Patented material is naturally adhering. Rinse clean to renew adhering properties.
- 3. Orbital-Genie For Standard Flask Clamps is equipped with a platform for accepting various size flasks or beakers. Clamps sold separately.



80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A.

Phone: (631) 567-4700 • Fax: (631) 567-5896

Toll Free: 888-850-6208

Customer Service: custsvc@scientificindustries.com

scientificindustries.com

Ordering Information and Accessories					
ORDERING INFORMATION					
Catalog No.	Description				
SI-1500	Orbital-Genie Shaker, Ratcheting Clamps, 120V				
SI-1600	Orbital-Genie Shaker, Standard Flask Clamps, 120V				
SI-1700	Orbital-Genie Shaker, Adhering Mat, 120V				
SI-1501	Orbital-Genie Shaker, Ratcheting Clamps, 230V, No Plug				
SI-1601	Orbital-Genie Shaker, Standard Flask Clamps, 230V, No Plug				
SI-1701	Orbital-Genie Shaker, Adhering Mat, 230V, No Plug				
SI-1502	Orbital-Genie Shaker, Ratcheting Clamps, 230V, European Plug				
SI-1602	Orbital-Genie Shaker, Standard Flask Clamps, 230V, European Plug				
SI-1702	Orbital-Genie Shaker, Adhering Mat, 230V, European Plug				
SI-1503	Orbital-Genie Shaker, Ratcheting Clamps, 230V, British Plug				
SI-1603	Orbital-Genie Shaker, Standard Flask Clamps, 230V, British Plug				
SI-1703	Orbital-Genie Shaker, Adhering Mat, 230V, British Plug				
SI-1504	Orbital-Genie Shaker, Ratcheting Clamps, 230V, Swiss Plug				
SI-1604	Orbital-Genie Shaker, Standard Flask Clamps, 230V, Swiss Plug				
SI-1704	Orbital-Genie Shaker, Adhering Mat, 230V, Swiss Plug				
SI-1506	Orbital-Genie Shaker, Ratcheting Clamps, 230V, Australian Plug				
SI-1606	Orbital-Genie Shaker, Standard Flask Clamps, 230V, Australian Plug				
SI-1706	Orbital-Genie Shaker, Adhering Mat, 230V, Australian Plug				
SI-1505	Orbital-Genie Shaker, Ratcheting Clamps, 100V				
SI-1605	Orbital-Genie Shaker, Standard Flask Clamps, 100V				
SI-1705	Orbital-Genie Shaker, Adhering Mat, 100V				



Flask Clamps (for SI-1600) are of stainless steel construction for added strength and safety. Mix and match different size clamps on same platform or all same size (capacity is indicated). Includes mounting screws. Note: Flask Clamps can also be used on SI-1700 by removing the blue Adhering

Double Sided Adhesive Tape and Adhering Mat (for SI-1600 or SI-1700) are alternatives to flask clamps for securing vessels to shaker platform. Both are ideal for containers such as bottles and dishes that cannot be held with conventional clamps. Note: Tape and Mat cannot be used on SI-1500.

ORDERING INFORMATION						
Catalog No.	Description	Platform Capacity	Catalog No.	Description		
SI-1610	100 mL Flask Clamp	36	SI-1616	Double Sided Adhesive Tape Roll, 60in x 1in (1.5m x 25.4mm)		
SI-1611	125/150 mL Flask Clamp	36	SI-1617	Adhering Mat, Blue, 18in x 18in (460mm x 460mm)		
SI-1612	250 mL Flask Clamp	25	SI-1510	Universal Ratcheting Clamp		
SI-1613 500 mL Flask Clamp		16				
SI-1614	1 L Flask Clamp 9			CLAUP - 2		
SI-1615	2 L Flask Clamp	5				

## **Specifications**

Speed10 - 300 RPM including stir reversePlatform Dimensions (D x W)460 x 460mm (18 x 18in)Timer1 min - 99 hours w/Alarm or ContinuousOverall Dimensions (D x W)590 x 500mm (23 x 20in)Orbit35mmMaximum Load12 Kg (26 lbs.)Selectable Ramp Rates10 - 100 rpm/secWeight21 Kg (46 lbs.)







