



**BIOHIT**

Innovating for Health

***Limited time offer!***

**Purchase any *e*LINE<sup>®</sup>  
and receive a Charging Stand and  
3 Ergomate microplate holders FREE**


***Over a  
\$230 value FREE!***



***Valid September 1, 2010 through December 31, 2010  
See back for additional information and ordering details***

# BIOHIT

Innovating for Health



## eLINE Performance specifications and ordering details

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Vol. (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Filters		Tip ID
							Standard	Plus	
730020	1-channel	0.2-10 µl	0.05	10 5 1 0.2	0.90 % 1.00 % 2.50 % 12.00 %	0.40 % 0.70 % 1.50 % 10.00 %	-	-	A, L, M
730040	1-channel	5-120 µl	0.50	120 60 12 5	0.40 % 0.60 % 2.00 % 3.50 %	0.15 % 0.20 % 1.00 % 1.50 %	721008	721018	C, D, V
730060	1-channel	10-300 µl	1.00	300 150 30 10	0.40 % 0.60 % 1.50 % 3.50 %	0.15 % 0.20 % 0.80 % 2.00 %	721007	721017	D, Q*
730080	1-channel	50-1000 µl	5.00	1000 500 100 50	0.40 % 0.60 % 1.50 % 2.00 %	0.15 % 0.20 % 0.50 % 1.00 %	721006	721016	E, G, S
730100	1-channel	100-5000 µl*	10.0	5000 2500 500	0.50 % 0.80 % 1.00 %	0.15 % 0.20 % 0.40 %	721006	721016	J
730320	8-channel	0.2-10 µl	0.05	10 5 1	0.90 % 1.50 % 4.00 %	0.50 % 0.80 % 3.00 %	-	-	A, L, M
730340	8-channel	5-120 µl	0.50	120 60 12	0.50 % 0.70 % 2.00 %	0.20 % 0.30 % 1.50 %	721008	721018	C, D, V
730360	8-channel	10-300 µl	1.00	300 150 30	0.50 % 0.70 % 2.00 %	0.20 % 0.30 % 1.00 %	721007	721017	D, Q*
730390	8-channel	50-1200 µl	5.00	1200 600 120	0.50 % 1.00 % 2.50 %	0.20 % 0.30 % 1.00 %	721006	721016	H, T, Z
730420	12-channel	0.2-10 µl	0.05	10 5 1	0.90 % 1.50 % 4.00 %	0.50 % 0.80 % 3.00 %	-	-	A, L, M
730440	12-channel	5-120 µl	0.50	120 60 12	0.50 % 0.70 % 2.00 %	0.20 % 0.30 % 1.50 %	721008	721018	C, D, V
730460	12-channel	10-300 µl	1.00	300 150 30	0.50 % 0.70 % 2.00 %	0.20 % 0.30 % 1.00 %	721007	721017	D, Q
730490	12-channel	50-1200 µl	5.00	1200 600 120	0.50 % 1.00 % 2.50 %	0.20 % 0.30 % 1.00 %	721006	721016	H, T, Z

\* Note: Minimum volume in P-mode is 500 µl. The 100 µl is possible in d-mode. The specifications are type test specifications achieved under strictly controlled conditions (ISO 8653). Due to Biohit's continuous R&D, specifications may change without prior notice.

**PIPETTING ACADEMY™**

WATCH VIDEO ON

[www.biohit.com/liquidhandling](http://www.biohit.com/liquidhandling)

To place an order please contact a participating Biohit dealer or contact your Biohit representative by calling 1-800-922-0784 or email to: [pipet@biohit.com](mailto:pipet@biohit.com).

Please visit our website: [www.us.biohit.com](http://www.us.biohit.com) for more information. In order to receive your free goods, please visit: [us.biohit.com/redemption.html](http://us.biohit.com/redemption.html) to download your redemption form. Once you fill out the form, you may send it to Biohit via email or fax for processing. Please allow up to 4 weeks for delivery of free goods. Offer cannot be combined with any other special promotions. Multiple redemptions permitted. Promotion not valid where prohibited, and may be discontinued without notice.

