



Pharmacy Refrigerators

MBI684



Audible and Visual Alarms

Audible buzzer and visible
flashing light



Multiple Alarm Functions

High/low temperature,
power failure, door ajar,
sensor error, low battery
remote system.



Air Cooling, Auto Defrost

Better cabinet temperature
uniformity, auto defrost

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

MBI684

Technical Data

Cabinet Type	Upright,solid door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	HC
Noise((dB(A))	41

Performance

Temp Range(°C)	2~8
-----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220~240/50	115/60
Power(W)	84	150
Electrical Current (A)	0.65	1.7

Dimensions

Capacity(L/Cu. Ft)	68/2.4
Net/Gross Weight(approx)	38/40(kg) 83.8/88.2(lbs)
Interior Dimensions(W*D*H)	415*385*505(mm) 16.3*15.2*19.9(in)
Exterior Dimensions(W*D*H)	495*580*660(mm) 19.5*22.8*26.0(in)
Packing Dimensions(W*D*H)	540*640*710(mm) 21.3*25.2*28.0(in)
Container Load(20'/40'/40'H)	108/228/228



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves /Basket	2/1

Others

Certificate	CE	UL
-------------	----	----





Pharmacy Refrigerators

MBI684



6

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Audible and Visual Alarms

Audible buzzer and visible
flashing light



Multiple Alarm Functions

High/low temperature,
power failure, door ajar,
sensor error, low battery
remote system.



Air Cooling, Auto Defrost

Better cabinet temperature
uniformity, auto defrost





Pharmacy Refrigerators

MBI684



6

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Audible and Visual Alarms

Audible buzzer and visible
flashing light



Multiple Alarm Functions

High/low temperature,
power failure, door ajar,
sensor error, low battery
remote system.



Air Cooling, Auto Defrost

Better cabinet temperature
uniformity, auto defrost

