




Pharmacy Refrigerators MBI9404



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Heated Double Pane Glass Door

Clear display of stored products without condensation
at 32°C and 85% relative humidity condition



Self-Closing Door

Energy conserving · maintains temperature
uniformity to protect valuable samples

Specifications

MBI9404

Technical Data

Cabinet Type	Upright, double glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	45

Performance

Temp Range(°C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

115/60

Power Supply(V/Hz)	220-240/50	208~230/60	
Power(W)	850	770	850
Electrical Current(A)	4.5	4.5	11

Dimensions

Capacity(L/Cu.Ft)	890/31.4
Net/Gross Weight(approx)	227/279(kg) 500.4/615.1(lbs)
Interior Dimensions(W*D*H)	1030*590*1425(mm) 40.6*23.2*56.1(in)
Exterior Dimensions(W*D*H)	1130*755*1980(mm) 44.5*29.7*78.0(in)
Packing Dimensions(W*D*H)	1250*850*2130(mm) 49.2*33.5*83.9(in)
Container Load(20'/40'/40'H)	8/18/18



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	12/-
USB Interface	Optional

Others

Certificate	CE
-------------	----