



Smart alarm and alerts for abnormal temperature, low storage, deadline for clinical trials

Net Weight / Gross Weight (Kg)	Shelf
116/139	5

Intelligent Management Plan byHaier U'Reagent


Cloud-shared data Rapid and smart inventory system

Instant generation of summary data report, clear transparency of transaction

RFID based intellectual technology tracks real time information on stored products. Effective management of medicines and reagents allows information of storage inventory instantly available at a stroke of the keyboard. Cloud based big data supports shared platform, permitting manufacturers, hospitals, quality control exercise real time monitoring. Technical data thus becomes the origin of productivity.



and smartly controlled. (It is suitable for hospitals and third party testing organizations)



Manpower savings

Smart process control improves accuracy and efficiency, saving manpower for reagent management.



Optimized procurement process reduces waste

Control reagent usage timely. Procure supply smartly to reduce shortage and reduce inventory.



Rigorous reagent management to reduce risk

Constant temperature for reagent storage
Automatically trigger temperature alarms when temperature reaches critical temperature, ensuring safety of reagent usage.



Big data platform for information clarity



Authorized usage for effective management

Set up different levels of authority based on responsibility



Remote control of terminal equipment with visible operation status

Real time inquiry of all terminal units

