



Pharmacy Refrigerator

MBI390F4



Frequency Conversion Energy Saving

Variable frequency compressor with optimized cooling system, more energy efficient.



Environmentally Friendly

Hydrocarbon compressor and refrigerant, the refrigerant does not contain any fluorine or chlorine which is super friendly to the environment.



Large Space, Superior Uniformity

The space is efficiently used, and the utilization rate is high and the temperature uniformity is $\pm 3^{\circ}\text{C}$.



Low Noise

Optimized system and noise reduction design to reduce the noise level to dB(A) (220V/50Hz).



Stable performance, multiple alarm modes

Three alarm modes (light, sound, remote); six types of alarm.



Microprocessor Control System



Self Closing Door with 90° Stay Open Feature.



Precise Temperature Control



Security lock



Optional USB and Printer



Scope of Application



This product is suitable for pharmacy, pharmaceutical manufacturers, anti-epidemic stations, health centers and hospitals as key examples, which is used to store biological products and articles that need to be stored at 2~8°C.

Specifications



MBI390F4

Technical Data		
Cabinet Type	Upright, Single door	
Climate Class	N	
Cooling Type	Forced air cooling	
Defrost Mode	Auto	
Refrigerant	HC	
Noise((dB(A))	43	43

Performance	
Temp Range(°C)	2~8

Control	
Controller	Microprocessor
Display	LED

Electrical Data		
Power Supply(V/Hz)	220~240/50&60	115/60
Power(W)	210	115
Electrical Current(A)	1	1.5

Dimensions		
Capacity(L/Cu.Ft)	106/129	
Net/Gross Weight(approx)	(kg)	255.8/306.6
	(lbs)	530*555*1380
Interior Dimension(W*D*H)	(mm)	20.9*21.9*54.3
	(in)	665*710*1965
Exterior Dimension(W*D*H)	(mm)	26.2*28.0*77.4
	(in)	690*790*2110
Packing Dimension(W*D*H)	(mm)	27.1*31.0*82.9
	(in)	21/42/42
Container Load(20'/40'/40'H)		



Accessories		
Caster	Y	
Foot	Y	
Porthole	Y	
Shelves/Drawers	7/-	
USB Interface	Optional	Y
Temperature Recorder	Optional	

Alarms	
High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Other		
Certificate	CE	UL