

65mm RODAC Dishes

Features

1. High clarity polystyrene
2. Transparent surface
3. Numeric mark on the bottom for easy identification
4. PA/PE packaging
5. Vent for gas exchange
6. Stackable for easy storage and handling
7. DNase/RNase and Human DNA Free, Non-pyrogenic
8. Sterilized by E-Beam, SAL=10⁻⁶

Description

Temperature Range Approx.-86 to+64℃
Shelf Life : 3 Years after date of production (ensure package is in good)
Sterile : Yes



Product Range

Cat.No.	Spec (mm)	Volume (mL)	/Pack	/Case
722011	65*15	13-15	20	500



Technical Drawing

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

Cap



Bottom



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA