

MBI SCI280-Pro

LED Digital Magnetic Hotplate Stirrer

Features:

- Digital control system is ideal for applications of Analytical and biological Laboratory.
- Temperature to 280°C
- Electronic speed control, the max speed is 1500rpm.
- Stainless steel with ceramic coated hotplate allows immediate heat transfer for quick heat-up times, chemical-resistant and scratchproof
- External PT1000 temperature sensor accessory is available with control accuracy of $\pm 0.5^{\circ}\text{C}$.
- LED display shows temperature and speed.
- 50°C HOT warning indicates residual hotplate temperature.
- ABS casing with heat resistant and resistance of weak acids and alkalis



MBI SCI280-Pro
temperature sensor
PT1000 accessory

