



ORBi-SHAKER™ XL

MAXIMUM CAPACITY, BENCHTOP MODEL



CE



with **MAGiC CLAMP™**
PATENTED

- + **High Capacity, benchtop model**
- + **Circular 19mm orbit for aeration and mixing**
- + **Exchange flask clamps in one second with MAGic Clamp™ platform.**

Technical Data

Speed Range:	30 to 300 rpm
Speed Increment:	1 rpm
Orbit:	19mm (3/4")
Max. Capacity:	5 x 2L or 9 x 1L Flasks
Max. Load:	65 lbs (30kg)
Platform Dimensions:	18 x 18 in. / 45 x 45 cm (Included)
Timer:	1 min. to 99 hours / continuous
Weight:	65 lbs / 30 kg
Dimensions:	18.5 x 21.25 x 5 in. / 47 x 54 x 12.5 cm
Electrical:	120V or 230V, 50-60Hz

The Orbi-Shaker XL offers flask capacity that far exceeds most benchtop models. With a large 18 x 18 inch platform and a powerful, beltless motor drive, the Orbi-Shaker XL is capable of mixing up to five 2L flasks or nine 1L flasks simultaneously. The included, non-slip rubber mat platform is ideal for low speed applications, while the convenient MAGic Clamp™ platform is also available for high speed mixing of Erlenmeyer flasks, tube racks or microplates. The MAGic Clamp™ allows for instant flask clamp exchange without the use of tools, screws or any other hardware.

Ordering Information

Orbi-Shaker XL, High capacity, benchtop shaker
Includes 18x18" platform with non-slip mat, 115V

Optional Platforms and Accessories:

MAGic Clamp universal platform, (18x18")
 Flask Clamp, 50ml (max. 49)
 Flask Clamp, 125ml (max. 36)
 Flask Clamp, 250ml (max. 25)
 Flask Clamp, 500ml (max. 16)
 Flask Clamp, 1000ml (max. 9)
 Flask Clamp, 2000ml (max. 5)
 Microplate Clamp (max. 8)
 Tube Rack, 20x50 & 30x15 ml tubes (max. 4)
 Tube Rack, 32x15 tubes (max. 8)
 Tube Rack, 15x50 tubes (max. 8)
 Bottle Rack, 4x500/600ml conical bottles (max. 4)
 Tube Rack, 12 x 1.5/2.0ml horizontal (max. 6)
 Tube Rack, 4 x 15ml horizontal (max.8)
 Tube Rack, 2 x 50ml horizontal (max.8)
 Adjustable flask clamp for MAGic clamp platform, accepts 125ml (max. 36), 250ml (max. 25) or 500ml (max. 16)
 Spring loaded platform, for uncommon vessels
 6 x microplate platform (standard or deepwell)

*To order item in 230V, add (-E) to item number.

