

MBI-spense Bottle-top Dispensers



Quick and easy volume adjustment

Smooth piston design provides effortless dispensing

Closure cap protects
against reagent contact and is easy to attach
and remove, even while wearing gloves.



Adapters allow fitting to standard
reagent bottles

Ordering Information

Cat No.	Volume	Increments	Accuracy	Precision
MBI62009990	0.5-5ml	0.10ml	0.5%	0.1%
MBI62009991	1.0-10ml	0.20ml	0.5%	0.1%
MBI62009992	2.5-25ml	0.50ml	0.5%	0.1%
MBI62009993	5.0-50ml	1.00ml	0.5%	0.1%

Accessories

Cat No.	Description
MBI06399926	1Ltr Autoclavable Amber Bottle (45mm neck)
MBI06399926	2.5Ltr Autoclavable Amber Bottle (45mm neck)
MBI06399958	S40 Spare Adapter (40mm)
MBI06399906	GL32 Spare Adapter (32mm)
MBI06399907	GL38 Spare Adapter (38mm)
MBI06399908	GL25 Spare Adapter (25mm)
MBI06399919	GL28 Spare Adapter (28mm)

Features:

- Excellent chemical resistance
- Entirely autoclavable at 121°C
- Four ranges of bottle top dispenser covers a volume range from 0.5ml to 50ml
- Easy to clean and maintain
- Made of PTFE, ETFE, FEP, BSG, PP
- 2.2g/cm³ vapor pressure Max. 500mbar, viscosity max. 500mm²/s, temperature max. 40°C, density max. 2.2g/cm³
- Supplied with S40, GL32, GL38, GL25, GL28 sized adapters.
- **3 YEAR WARRANTY**