

# Complete solutions for your set-up from a single supplier

## Set-up solutions for one rotary evaporator

SCOPE OF DELIVERY\*

Rotary Evaporator

+ Vacuum Controller

+ Vacuum Pump

+ Chiller



RV 3 pro V Complete



RV 8 pro V Complete



RV 10 digital pro V Complete



RV 10 auto pro V-C Complete

	RV 3 pro V Complete	RV 8 pro V Complete	RV 10 digital pro V Complete	RV 10 auto pro V-C Complete
Ident. No.	0010011410	0010011963	0010011411	0010007353
Manual set-up	✓	✓	✓	✓
Heating bath temperature of 180 °C	–	✓	✓	✓
Motorized lift	–	–	✓	✓
Automatic boiling point detection and controlled lift/heating bath and cooling water management	–	–	–	✓

## Set-up solution for two rotary evaporators

SCOPE OF DELIVERY

> 2 x RV 10 digital rotary evaporator

+ 1 x VACSTAR digital vacuum pump

+ 1 x VSE 1 vacuum safety emission condensor

+ 2 x VC 10 pro vacuum controller

+ 1 x RC 2 GREEN basic chiller

+ 2 x RV 10.5005 set of hoses

+ 2 x VC 10.300 check valve

+ 2 x VC 10.400 table stand

+ 2 x Y-hose connector, 8 mm



\* The Complete Solutions include all required tubes and connection cables.  
Details on individual components can be found on the website.