

Cluster tube

NEST 1.1mL 96-well Cluster Tube Box conforms to the ANSI-SBS standard and is compatible with multi-channel pipettes and automated liquid workstations. The raw material of tube body is imported PP, which meets the USP-CLASS 6 standard. It has high transparency and good chemical stability, and is suitable for low-temperature storage of cells, drugs, tissues, sera, and other reagents.

Complete Specifications





- NEST 1.1mL Cluster Tubes are available in single tube, 8-strip, and 12-strip, and the tube caps are available in 8-strip and 12-strip for customers to choose from.
- They are compatible with multi-channel pipettes and automated liquid workstations with SBS standard size carriers.

Exquisite Design

- The one-way cap design, with letter and number grids, ensures positive identification of samples.
- The bottom and side walls of the box holder are hollowed out for good thermal conductivity, and the box supports direct immersion in a water bath.

High Quality Raw Materials

- The raw material is imported PP, which conforms to USP class 6 standards. The 8-strip and 12-strip caps are made of TPE, which has better sealing performance.
- Resistant to low temperature storage at -80°C, and autoclavable.




Cluster Tube

Schematic Diagram	Name	Packaging			Material	Cat.NO. (sterile)	Cat.NO. (non-sterile)
	1.1mL Single Cluster Tube	Bagged	960 pcs/bag	5 bags/case	PP	/	628111
		Boxed	96 pcs/box 10 boxes/pack	5 bags/case	PP	628114	628112
	1.1mL 8-strip Cluster Tube	Bagged	120 pcs/bag	5 bags/case	PP	/	628121
		Boxed	12 pcs/box 10 boxes/pack	5 bags/case	PP	628124	628122
	1.1mL 12-strip Cluster Tube	Bagged	80 pcs/bag	5 bags/case	PP	/	628131
		Boxed	8 pcs/box 10 boxes/pack	5 bags/case	PP	628134	628132


Cluster Tube Cap

Schematic Diagram	Name	Packaging			Material	Cat.NO. (sterile)	Cat.NO. (non-sterile)
	8-strip Cluster Tube Cap	Bagged	120/bag	5 bags/case	TPE	628923	628921
	12-strip Cluster Tube Cap	Bagged	80/bag	5 bags/case	TPE	628933	628931


● Technical Drawing




1.1mL Single Cluster Tube




1.1mL 8-strip Cluster Tube






1.1mL 12-strip Cluster Tube



8-strip Cluster Tube Cap



12-strip Cluster Tube Cap



Technical Drawing

Box Cover