

15mL Centrifuge Tubes

Features

1. Made from virgin polypropylene to meet the general requirements of the laboratory.
2. Printed clear graduations with marking area
3. PP for the tube and HDPE for the cap
4. One-hand operation design
5. Sterilized by E-Beam
6. DNase, RNase and pyrogen free


Description

Temperature Range : Approx.-80 to+121℃
 Shelf Life : 3 years after date of production
 Sterile : Yes

Product Range

Cat.No.	Volume (mL)	Packaging	Max.RCF (xg)	/Pack	/Case
601051	15	PS Rack	12000	50	500
601052		Bagged			

Technical Drawing



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA