

50mL Serological Pipettes

NEST serological pipettes can be used in both tissue culture and bacterial culture experiments.

Features


1. Made from high clarity polystyrene.
2. Individually wrapped
3. Edged color-coded band for identification
4. Negative graduations for extra volume
5. Sterilized by E-beam
6. Non-pyrogenic, non-cytotoxic and non-haemolytic
7. DNase/RNase free
8. Polyolefin fiber filter to reduce contamination

Description

Temperature Range : Approx.0 to+40℃
 Shelf Life : 4 years after date of production
 Sterile : Yes
 Capacity tolerance : $\pm 2\%$

Product Range

Cat.No.	Volume (mL)	Length (mm)	/Pack	/Box	/Case
329001	50	344	1	100	600



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA