

402301 PCR Plates

For PCR machine whose fluorescent signals is read at the top, Nest white PCR plates can effectively prevent the cross interference between different holes. In addition, Nest white PCR plate can read low fluorescent signals with high stability. Nest transparent PCR plates apply to the PCR machines whose fluorescent signals is read at the bottom.

Features

1. Made from high quality polypropylene to ensure the minimal loss of reaction solution.
2. Flat surface, thick and solid, not easy to deform
3. The black printing mark is designed for easy identification
4. Compatible with NEST pressure sensitive film, self-adhesive film and hot sealing film
5. Autoclavable

Description


Temperature Range : Approx.-80℃ to+121℃
Shelf Life : 5 Years after date of production (ensure package is in good)
Sterile : No

Product Range

Cat.No.	Volume (mL)	Skirt Type	Cut Position	Color	/Pack	/Case
402301	0.2	Semi Skirt	A12	Transparent	25	100

Technical Drawing

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA