

# PCR Plates

For PCR machine whose fluorescent signals is read at the top, Nest white PCR plates can effectively prevent the cross interference between different holes. In addition, Nest white PCR plate can read low fluorescent signals with high stability. Nest transparent PCR plates apply to the PCR machines whose fluorescent signals is read at the bottom.

## Features

1. Made from high quality polypropylene to ensure the minimal loss of reaction solution.
2. Flat surface, thick and solid, not easy to deform
3. The black printing mark is designed for easy identification
4. Compatible with NEST pressure sensitive film, self-adhesive film and hot sealing film
5. Autoclavable

## Description

Temperature Range : Approx.-80°C to+121°C  
Shelf Life : 5 Years after date of production ( ensure package is in good )  
Sterile : No

## Product Range

Cat.No.	Volume (mL)	Skirt Type	Cut Position	Color	/Pack	/Case
402011	0.2	No Skirt	H1	White	25	100

## Technical Drawing

<b>NEST China</b>	Tel	: 0510-6878 8698 / 6878 8718	<b>NEST USA</b>	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA