

ELISA Plates

Features

1. PS material, designed for ELISA method
2. Conform to SBS/ ANSI
3. Both detachable and undetectable plates are available
4. Flat and transparent bottom
5. Excellent well uniformity, CV <5%
6. DNase/RNase and Human DNA Free, Non-pyrogenic


Description

Temperature Range : Approx.-86 to+64℃
 Shelf Life : 3 years after date of production
 Sterile : Yes

Product Range

Cat.No.	Description	Spec (well)	Color	/Pack	/Case
504201	Detachable, High Binding, Clear	12*8	White bottom & Transparent well	5	50

Technical Drawing


**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.230, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018