

0.5ml Cryogenic Vials

Features

1. USP VI Polypropylene
2. Non-pyrogenic, DNase/RNase and human DNA free
3. Round Bottom, Self-standing
4. Solid, one-piece construction ensures the barcode not fall off during the storage under low temperature
5. Autoclavable
6. Sterilized by E-beam, SAL=10⁻⁶
7. Clear graduations and white marking area
8. Recommend to store in gas phase of the nitrogen tank
9. Internal thread Cap: Special screw cap designed with silicone O-ring to secure the sealing under various condition.
10. External thread Cap: Special screw cap designed without silicone ring to avoid the contamination risk.


Description

Temperature Range	:	Approx.-196°C to 121°C
Shelf Life	:	5 Years after date of production (ensure package is in good)
Sterile	:	Yes

Product Range

Cat.No.	Volume (mL)	Cap Style	Packaging	/Pack	/Case
618001	0.5	External thread	Plastic bag	50	2,000
618101		Internal thread			
618003		External thread	9*9/Box	81	1,296
618103		Internal thread			
618102		Internal thread	10*10/Box	100	1,600

Technical Drawing


**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.230, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018