

3mL Pasteur Pipettes

It is suitable for drawing, transferring and blending a certain volume of liquid. Both sterile and non-sterile Pasteur pipettes are available.

Features

1. NEST Pasteur pipettes are made of LDPE and manufactured in ISO 8 grade cleaning room.
2. Good chemical stability
3. Transparent with clear graduations for easy observation
4. Pipette can be twisted for use with micro vessels and special-shaped containers
5. High precision and good repeatability
6. Individually packed pasteur pipettes are sterilized by EO.
7. DNase/RNase, endotoxin and pyrogenic free


Description

Temperature Range : Approx.0 to+40℃


Shelf Life : 2 years after date of production

Product Range

Cat.No.	Volume (mL)	Length (mm)	Pack	Sterile	/Pack	/Box	/Case
318212	3	162	individually packed	Yes	1	500	2,000
318314	3	182				200	2,000
318231	3	162	Bagged	No	/	500	5,000



318212/318231



318314

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA