

0.2mL PCR Tube

Features


1. Made from USP grade VI polypropylene
2. Compatible with most thermal cyclers and real-time PCR machine
3. Extra-thin walls provide the optimal thermal transfer and shorter cycle times
4. Easy open leak proof cap makes the loss of reaction less than 5%
5. Free of Dnase and Rnase, human DNA and PCR inhibitor

Description

Temperature Range : Approx.-80°C to+121°C
 Shelf Life : 5 Years after date of production (ensure package is in good)
 Sterile : No

Product Range

Cat.No.	Volume (mL)	Specification	Cap Type	Color	/Pack	/Case
401001	0.2	High Profile	Flat Cap	Transparent	1,000	10,000



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018