

PCR 8-Strip Tube Caps

Features


1. Made from USP grade VI polypropylene.
2. Compatible with both NEST PCR 8-strip tube and PCR plate.
3. Domed and flat caps are available
4. Tight sealing, easy to open
5. Flat cap is ideal for qPCR
6. The design of domed cap can better prevent liquid evaporation

Description

Temperature Range : Approx. -80°C to +121°C
 Shelf Life : 5 Years after date of production (ensure package is in good)
 Sterile : No

Product Range

Cat.No.	Cap Type	Color	/Pack	/Case
406012	Flat	Transparent	125	1250
406112	Domed	Transparent	125	1250



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA