

# 5mL Serological Pipettes

NEST serological pipettes can be used in both tissue culture and bacterial culture experiments.

## Features


1. Made from high clarity polystyrene.
2. Individually wrapped
3. Edged color-coded band for identification
4. Negative graduations for extra volume
5. Sterilized by E-beam
6. Non-pyrogenic, non-cytotoxic and non-haemolytic
7. DNase/RNase free
8. Polyolefin fiber filter to reduce contamination

## Description

Temperature Range : Approx.0 to+40℃  
 Shelf Life : 4 years after date of production  
 Sterile : Yes

## Product Range

Cat.No.	Volume (mL)	Length (mm)	/Pack	/Box	/Case
326001	5	340	1	200	800



5ml

<b>NEST China</b>	Tel	: 0510-6878 8698 / 6878 8718	<b>NEST USA</b>	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018