



## **More Than Just Vortexing!**

## Another revolutionary product from the makers of the Vortex-Genie!

- True pulsing action provides more collisions and random motion.
- For applications requiring more aggressive action than standard vortexing.
- Variable "RUN/STOP" time provides quick pulsing or interval mixing.
- Fully programmable "ON" and "REST" function reduces heat buildup within the vessel, while keeping materials in suspension. Ideal for heat sensitive samples.
- O Digital control and display of speed, time, and "RUN/STOP" time, for reproducibility and accuracy, with true speed display.
- Use with all Vortex-Genie accessories.
- May also be used as a standard "continuous" Vortex Mixer.
- Oldeal for applications such as re-suspensions of lyophilized materials, creating solutions, enzymatic reactions and much more.

Ordering	Information
----------	-------------

Catalog No.	Description
SI-P236	Vortex-Genie Pulse, 120V
SI-P246	Vortex-Genie Pulse, 230V, No Plug
SI-P256	Vortex-Genie Pulse, 230V, Euro Plug
SI-P266	Vortex-Genie Pulse, 230V, British Plug
SI-P276	Vortex-Genie Pulse, 230V, Swiss Plug
SI-P296	Vortex-Genie Pulse, 240V, Australian Plug
SI-P286	Vortex-Genie Pulse, 100V

## **Specifications**

Supplied 3-inch Platform, Pop-off Cup
Speed 500 - 3000 RPM (2850 for 50Hz models)
Timer - Hands Free 1-99 minutes or continuous
Pulse - ON/OFF 0.1 seconds - 99 minutes
Weight 4 kgs (8.8 lbs)
Dimensions (D x W x H) 165 x 122 x 165mm

(6.5 x 4.8 x 6.5in)









80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A.

Phone: (631) 567-4700 • Fax: (631) 567-5896

Toll Free: 888-850-6208

Customer Service: custsvc@scientificindustries.com

scientificindustries.com

## Multi-Task Genie™ Family Accessories

Accessories for all Vortex-Genie® 2 mixers and Vortex-Genie® Pulse



Fig. 1 Model H301 Multiple Sample Starter Set





Fig. 4 Model H302 Large Sample Set

Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
1	504-0233-00	Model H301 Multiple Sample Starter Set	13	SI-0564	TurboMix™ Attachment
	146-6005-00	6-inch Platform	14	SI-0565	Universal Microtube Holder*
2	504-0234-00	Microtube Foam Insert (2)**	15	SI-0511	Large Ampule/Tube Attachment
3	504-0235-00	Microplate Foam Insert (2)**	16	SI-0570	Small Ampule/Tube Attachment
3A	504-0232-00	Deepwell Microplate Foam Insert (2)**	17	SI-H506	Horizontal 50ml Tube Holder (6 tubes)
4	503-0282-00	Model H302 Large Sample Set	18	SI-H512	Horizontal 15ml Tube Holder (12 tubes)
5	504-0039-00	Recessed Platform	19	SI-H524	Horizontal Microtube Holder (24 tubes)
6	503-0278-00	9-16mm Tube Foam Inserts (2)***	20	SI-V203	Vertical High Speed 50ml Tube Holder (3 tubes)
7	503-0279-00	14-29mm Tube Foam Inserts (2)***	21	SI-V206	Vertical High Speed 15ml Tube Holder (6 tubes)
8	503-0280-00	29-37mm Tube Foam Inserts (2)***	22	SI-V209	Vertical High Speed Microtube Holder (9 tubes)
9	503-0277-00	Blank Foam Inserts (2)***	23	SI-V506	Vertical 50ml Tube Holder (6 tubes)
	503-0283-00	Assorted Foam Inserts (4)***	24	SI-V512	Vertical 15ml Tube Holder (12 tubes)
10	568-0001-00	Elastic Bands (2)***	25	SI-V524	Vertical Microtube Holder (24 tubes)
11	146-3011-00	Pop-off Cup	26	SI-V525	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
12	0K-0500-902	3-inch Platform with Rubber Cover		0K-0236-408	Bottom Closure with Screws Kit †
	580-2013-00	Rubber Cover for 3-inch Platform			

<sup>\*</sup> May be used with TurboMix attachment or Disruptor Genie®

† May not fit older models. Please inquire. NOTE: Vortex-Genie 1 may only be used with Pop-off Cup and 3-inch Platform

U.S. Pat. 5,707,861



80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A.

Phone: (631) 567-4700 • Fax: (631) 567-5896

Toll Free: 888-850-6208

Customer Service: custsvc@scientificindustries.com

scientificindustries.com

<sup>\*\*\* 6-</sup>inch Platform Required
\*\*\* Recessed Platform Required