

# Power Supply











MBIMP-310 / MBIMP-320 / MBIMP-510



Stackable case



Four pairs of outlet terminals

	MBIMP-100	MBIMINI-300	MBIMINI-500	MBIMP- 10	MBIMP-	MBIMP- 10
Cat. No.		 CE ENEC	 CE ENEC	 FC CE ENEC	 FC CE	 FC CE
Display	N/A	LED		2.4" TFT LCD		
Output voltage / Inc.	50V/100V	10-300V/1V	10-500V/1V	5-300V/1V		5-500V/1V
Output current / Inc.	400mA	10-400mA/1mA	1-400mA/1mA	1-700mA/1mA	10-3000mA/10mA	1-800mA/10mA
Output power / Inc.	40W	60W	120W	Output: 150W Input: 200W	300W	
Output type	Constant voltage	Constant voltage/Current		Constant voltage/Current /Power		
Timer	N/A	999 (min) with alarm/Continuous		Constant mode:9999 (min) with alarm/Continuous Programmable mode: 999 (min) with alarm/Continuous		
Rated voltage	100V-120V~ Only	100V-240V~				
Program	N/A	2-steps (V, mA,  )		Pre setting: Up to 6-steps (V, mA, W,  ) , 30 programmed files		

