

Shakers

DLAB rocking and 3D shakers provide perfect solution in wide fields including electrophoresis gel dyeing/ discoloration, sample washing, molecular hybridization, immunoprecipitation, Western blot and cell culture.

MBI SK-R330-Pro LCD Digital Rocking Shaker

Features:

- Rocking motion angle 9°
- Max. load capacity is 10kg with platform
- Speed range of 10~70rpm
- Double LCDs for independent display time and speed
- Various of platforms, accessories



MBI SK-R1807-S LED Digital Rocking shaker

Features:

- Rocking motion angle 7°
- Max. load capacity is 3kg with platform
- Speed range of 10~80rpm
- Double LEDs for independent display time and speed
- Various of platforms, accessories



MBI SK-R1807-E Economical Rocking Shaker

Features:

- Rocking motion angle 7°
- Max. load capacity is 3kg with platform
- Speed range of 0~80rpm
- Various of platforms, accessories



Shakers



Timer



Wide range of accessories



Specifications

	MBI SK-R330-Pro	MBI SK-R1807-S	MBI SK-R1807-E
Tilt angle	9°	7°	7°
Max. load capacity (with platform)	10kg	3kg	3kg
Motor type	DC motor	DC motor	DC motor
Motor input power	40W	16W	16W
Motor output power	24W	10W	10W
Speed range	10~70rpm	10~80rpm	0~80rpm
Speed display	LCD	LED	-
Timer display	LCD	LED	-
Time settings range	1min~19h59min	1min~19h59min	-
Voltage	100~240V,50/60Hz	100~240V,50/60Hz	100~240V,50/60Hz
Power	50W	20W	20W
Weight	9kg	3.2kg	3.2kg
Permissible ambient temperature and humidity	5~40°C,80%	5~40°C,80%	5~40°C,80%
Dimension[WxDxH]	360x410x200mm	270x330x130mm	270x330x130mm
Protection class	IP21	IP21	IP21

