

Centrifuge Tubes

Features

1. USP VI Materials
2. Meet the requirements of analytical level
3. High strength non-toxic materials
4. Printed clear graduations with marking area
5. PP for the tube and HDPE for the cap
6. One-hand operation design
7. Sterilized by E-Beam
8. DNase, RNase and pyrogen free

Description

Temperature Range : Approx.-80 to+121℃
 Shelf Life : 3 years after date of production
 Sterile : Yes

Product Range

Centrifuge Tubes, USP VI Materials

Cat. No.	Volume(Ml)	Color	Packaging	Max.RCF(xg)	/Pack	/Case
601001	15	transparency	PS Rack	12000	50	500
601002	15	transparency	Bagged	12000	50	500
602001	50	transparency	PS Rack	12000	25	500
602002	50	transparency	Bagged	12000	25	500
601201	15	brown	PS Rack	12000	50	500
602201	50	brown	PS Rack	12000	25	500

Centrifuge Tubes

Cat. No.	Volume(mL)	Packaging	Max.RCF(xg)	/Pack	/Case
601051	15	PS Rack	12000	50	500
601052	15	Bagged	12000	50	500
602051	50	PS Rack	12000	25	500

**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**


Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. © Wuxi NEST Biotechnology Co., Ltd.


602052	50	Bagged	12000	25	500
602072	50 (Self-Standing)	Bagged	5000	25	500

Technical Drawing


15mL Centrifuge tubes/ USP VI Materials



50mL Centrifuge tubes/ USP VI Materials



50mL Self-Standing Centrifuge tubes/ USP VI Materials



**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA