

384 Well Cell Culture Plate

Features

1. High clarity, 100% virgin polystyrene
2. Sterilized by E-beam
3. Vacuum Plasma tissue culture treatment
4. Clear lot number for batch traceability
5. Markings of well coordinates available for 12, 24, 48 96 and 384 well plates
6. Non-Pyrogenic, DNase/RNase Free
7. The package box of Nest cell culture plates can be used as reservoirs

Description

Temperature Range : Approx.-86 to+64℃
 Shelf Life : 3 Years after date of production (ensure package is in good)
 Sterile : Yes

Product Range

Cat.No.	Spec (Well)	Cell Growth Area (cm ²)	TC	Rcommended Medium Volume (mL)	Color	packaging	/Pack	/Case
761001	384	0.11	TC	0.025-0.05	Transparent	Plastic box	1	100
761301					Black			
761601					White			
761011			Transparent					
761311			Black					
761611			White					
761002			TC		Transparent	plastic bag	10	100
761012			Non-Treated					

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA


NEST TÜVRheinland



L21-231111004
06140

NEST

Cell Culture Plate
Flat Bottom with Low Evaporation Lid
Tissue Culture are Treated

[illegible]

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018