

96 Well Cell Culture Plate

Features

1. High clarity, 100% virgin polystyrene
2. Sterilized by E-beam
3. Vacuum Plasma tissue culture treatment
4. Clear lot number for batch traceability
5. Markings of well coordinates available for 12, 24, 48 96 and 384 well plates
6. Non-Pyrogenic, DNase/RNase Free
7. The package box of Nest cell culture plates can be used as reservoirs

Description

Temperature Range : Approx.-86 to+64℃
 Shelf Life : 3 Years after date of production (ensure package is in good)
 Sterile : Yes

Product Range

Cat.No.	Spec (Well)	Cell Growth Area (cm ²)	TC	Rcommended Medium Volume (mL)	Bottom	packaging	/Pack	/Case
701001	96	0.32	TC	0.1-0.2	Flat	Plastic box	1	100
701101					U			
701201					V			
701011			Non-Treated	0.1-0.2	Flat			
701111					U			
701211					V			
701002			TC	0.1-0.2	Flat	plastic bag	10	100
701012			Non-Treated		Flat			

**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.230, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**


Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA



Excellent bottom uniformity



Special lid design ensures low evaporation





Frosted and noise-free bottom design for easy grip



Lot No. for traceability

Flat Bottom





**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
Email : info@nest-wuxi.com
Website : www.cell-nest.com
Add : No.230, xida road, new district, Wuxi, Jiangsu, China



**NEST
USA**

Tel : +1-732-381-0268
Email : info@nestscientificusa.com
Website : www.nestscientificusa.com
Add : 1592 Hart St., Rahway, NJ07065, USA

U Bottom



V Bottom



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA