

100mm Cell Culture Dishes

NEST cell culture dishes are ideal for all types of culturing where microscopic observation is required. Numeric marks on the bottom of each dish allow users to identify the location of cells.

Features

1. High clarity polystyrene
2. Flat transparent surface for distortion-free observation
3. Vacuum plasma TC treatment, excellent cell adherence
4. Stackable for easy storage and handling
5. Sterilized by E-Beam
6. Non-pyrogenic

Description



Temperature Range Approx.-86 to+64℃
 Shelf Life : 3 Years after date of production (ensure package is in good)
 Sterile : Yes

Product Range



Cat.No.	Spec (mm)	Cell Growth Area (cm ²)	Rcommended Medium Volume (mL)	/Pack	/Case	Remark
704004	100*17.7	55	11-16.5	5	300	/
704002				10	300	/
704001				20	300	/
704202				10	300	With Griping Ring
704201				20	300	With Griping Ring

Technical Drawing

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA



With Gripping Ring



**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
Email : info@nest-wuxi.com
Website : www.cell-nest.com
Add : No.230, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
Email : info@nestscientificusa.com
Website : www.nestscientificusa.com
Add : 1592 Hart St., Rahway, NJ07065, USA