

T25 Cell Culture Flask

NEST cell culture flasks are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum free or serum-containing cultures.

Features

1. High clarity medical grade polystyrene material
2. Sterilized by E-Beam SAL=10⁻⁶
3. Nonpyrogenic, DNase free, RNase free, Non-Cytotoxicity
4. Packaged in sterile, zip-sealable bags
5. Two types of cap style : vent caps & plug seal cap
6. Vent caps with 0.22um hydrophobic filters for gas exchange without contamination
7. Fostered writing and clear graduation
8. Clear lot number for batch traceability
9. Stackable

Description

Temperature Range Approx.-86 to+64℃
 Shelf Life : 3 Years after date of production (ensure package is in good)
 Sterile : Yes

Product Range


Cat.No.	Cap Style	TC Treated	Cell Growth Area	Volume	Recommended Medium Volume	/Pack	/Case
707001	Plug Seal Cap	Yes	25 cm ²	50 ml	5-7.5 ml	10	200
707003	Vent Cap	Yes					
707011	Plug Seal Cap	No					
707013	Vent Cap	No					

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018

Technical Drawing



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA