

5-layer Cell Culture Flask

NEST cell culture flasks are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum free or serum-containing cultures.

Features

1. High clarity medical grade polystyrene material
2. Sterilized by E-Beam SAL=10⁻⁶
3. Nonpyrogenic
4. Packaged in sterile bag
5. Two types of cap style : vent caps & plug seal cap
6. Vent caps with 0.22um hydrophobic filters for gas exchange without contamination
7. TC treated , good for cell attachment and growth.

Description

Cat.No.	: 731001、731002
Growth area	: 870 cm ²
Cap Style	: Plug Seal Cap (731001)、Vent Cap (731002)
Temperature Range	Approx.-86 to+64℃
Shelf Life	: 3 Years after date of production (ensure package is in good)
Sterile	: Yes






Product Range

Cat.No.	Growth area	Cap Style	Size (mm)			TC Treated	/Pack	/Case
			Length	Width	Height			
731001	870 cm ²	Plug Seal Cap	203.6	120.5	84.3	Yes	1	8
731002		Vent Cap						

Recommendations

Please be gently when use the 5-layer cell culture flask to avoid air bubbles. Air bubbles will lead the culture medium flow from the upper layer to the lower layer.

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA