

2mL Serological Pipettes

NEST serological pipettes can be used in both tissue culture and bacterial culture experiments.

Features


1. Made from high clarity polystyrene.
2. Individually wrapped
3. Edged color-coded band for identification
4. Negative graduations for extra volume
5. Sterilized by E-beam
6. Non-pyrogenic, non-cytotoxic and non-haemolytic
7. DNase/RNase free
8. Polyolefin fiber filter to reduce contamination

Description

Temperature Range : Approx.0 to +40°C
 Shelf Life : 4 years after date of production
 Sterile : Yes

Product Range

Cat.No.	Volume (mL)	Length (mm)	/Pack	/Box	/Case
325001	2	280	1	400	2,400



2ml

NEST China	Tel :	0510-6878 8698 / 6878 8718	NEST USA	Tel :	+1-732-381-0268
	Email :	info@nest-wuxi.com		Email :	info@nestscientificusa.com
	Website :	www.cell-nest.com		Website :	www.nestscientificusa.com
	Add :	No.230, xida road, new district, Wuxi, Jiangsu, China		Add :	1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018