

15mL Centrifuge Tubes

Features

1. USP VI Materials
2. Meet the requirements of analytical level
3. High strength non-toxic materials
4. Printed clear graduations with marking area
5. PP for the tube and HDPE for the cap
6. One-hand operation design
7. Sterilized by E-Beam
8. DNase, RNase and pyrogen free


Description

Temperature Range : Approx.-80 to+121℃
 Shelf Life : 3 years after date of production
 Sterile : Yes

Product Range

Cat.No.	Volume (mL)	Packaging	Max.RCF (xg)	/Pack	/Case
601001	15	PS Rack	12000	50	500
601002		Bagged			

Technical Drawing



NEST China	Tel :	0510-6878 8698 / 6878 8718	NEST USA	Tel :	+1-732-381-0268
	Email :	info@nest-wuxi.com		Email :	info@nestscientificusa.com
	Website :	www.cell-nest.com		Website :	www.nestscientificusa.com
	Add :	No.230, xida road, new district, Wuxi, Jiangsu, China		Add :	1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018