

1.5mL Micro Centrifuge Tube

Features



1. USP VI grade polypropylene material, free from heavy metal
2. RCF: 30000 xg
3. Special safelock cap design provides secure lid sealing
4. Clear graduations and frosted marking area
5. Easy opening and closing
6. Autoclavable
7. DNase/RNase free, endotoxin level less than 0.1EU
8. Both sterile and non-sterile tubes are available

Description

Temperature Range : Approx.-80°C to+121°C
 Autoclavability : Autoclave at 121°C/15psi
 615001 Shelf Life : 5 years after date of production
 615601 Shelf Life : 3 years after date of production

Product Range

Cat.No.	Volume (mL)	Bottom	Sterile	/Pack	/Case
615001	1.5	Conical	NO	500	5,000
615601			Yes	50	4,000



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA