

0.2mL PCR 8-Strip Tubes

NEST has both transparent and white PCR strip tubes, which respectively apply to the general PCR and qPCR reactions.

Features

1. Made from USP grade VI polypropylene
2. Low profile tubes are designed to minimize the condensate contamination, particularly for the overnight experiment, and also reduce the evaporation of the reaction solution.
3. Free of Dnase and Rnase, human DNA and PCR inhibitor
4. White strip tubes increase the fluorescent signal strength of real-time PCR.

Description


Temperature Range	: Approx.-80℃ to+121℃
Shelf Life	: 5 Years after date of production (ensure package is in good)
Sterile	: No
Maximum centrifugal force	: 4000g

Product Range

Cat.No.	Volume (mL)	Specification	Cap Type	Color	/Pack	/Case
403002	0.2	High Profile	No Cap	Transparent	125	1250
403012	0.2	High Profile	No Cap	White	125	1250
404001	0.2	High Profile	Flat Cap	Transparent	125	1250

Technical Drawing

403002 & 403012





NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018

404001



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA