

PCR Rack

Features

1. Made from USP grade VI polypropylene
2. Fit for fixing all NEST PCR products.
3. Stackable
4. Transparent polycarbonate lid for well identification
5. Blur Rack & Transparent Lid


Description

Temperature Range : Approx.-80℃ to+121℃
 Shelf Life : 5 Years after date of production (ensure package is in good)
 Sterile : No

Product Range

Cat.No.	Spec	Cap Type	/Pack	/Case
407001	8*12	Flat	1 Rack & 1 Lid	25 Rack & 25 Lid
407101	8*12	Flat	5 Racks & 1 Lid	25 Rack & 5 Lid

Technical Drawing


**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.230, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA