

90mm Petri Dishes

Features

1. High clarity polystyrene
2. Flat transparent surface
3. Numeric mark on the bottom for easy identification
4. PA/PE packaging
5. Vent for gas exchange
6. Stackable for easy storage and handling
7. DNase/RNase and Human DNA Free, Non-pyrogenic
8. Sterilized by E-Beam, SAL=10⁻⁶


Description

Temperature Range Approx.-86 to+64℃
 Shelf Life : 3 Years after date of production (ensure package is in good)
 Sterile : Yes

Product Range

Cat.No.	Spec (mm)	Volume (mL)	TC-Treated	/Pack	/Case
752031	90*15	10 x 4 compartments	No	20	500

Technical Drawing


NEST
China

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.230, xida road, new district, Wuxi, Jiangsu, China

NEST
USA

Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA