

514201 ELISA Plates

Features

1. PS material, designed for ELISA method
2. Conform to SBS/ ANSI
3. Both detachable and undetectable plates are available
4. Flat and transparent bottom
5. Excellent well uniformity, CV <5%


Description

Shelf Life : 3 years after date of production
Sterile : Non-sterile

Product Range

Cat.No.	Description	Spec (well)	Color	/Pack	/Case
514201	Undetachable, High Binding, Clear	12*8	Transparent	5	50

Technical Drawing



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018