

5.0ml Cryogenic Vials

Features

1. USP VI Polypropylene
2. Non-pyrogenic, DNase/RNase and human DNA free
3. Round Bottom, Self-standing
4. Solid, one-piece construction ensures the barcode not fall off during the storage under low temperature
5. Autoclavable
6. Sterilized by E-beam, SAL=10⁻⁶
7. Clear graduations and white marking area
8. Recommend to store in gas phase of the nitrogen tank
9. Internal thread Cap: Special screw cap designed with silicone O-ring to secure the sealing under various condition.
10. External thread Cap: Special screw cap designed without silicone ring to avoid the contamination risk.

Description

Temperature Range	:	Approx.-196°C to 121°C
Shelf Life	:	5 Years after date of production (ensure package is in good)
Sterile	:	Yes


Product Range

Cat.No.	Volume (mL)	Cap Style	Packaging	/Pack	/Case
609001	5.0	External thread	Plastic bag	50	1,000
609101		Internal thread			
609003		External thread	9*9/Box	81	486
609103		Internal thread			


NEST China	Tel	:	0510-6878 8698 / 6878 8718	NEST USA	Tel	:	+1-732-381-0268
	Email	:	info@nest-wuxi.com		Email	:	info@nestscientificusa.com
	Website	:	www.cell-nest.com		Website	:	www.nestscientificusa.com
	Add	:	No.230, xida road, new district, Wuxi, Jiangsu, China		Add	:	1592 Hart St., Rahway, NJ07065, USA

Technical Drawing

External thread



Internal thread



**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
Email : info@nest-wuxi.com
Website : www.cell-nest.com
Add : No.230, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
Email : info@nestscientificusa.com
Website : www.nestscientificusa.com
Add : 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018