

1mL Pasteur Pipettes

It is suitable for drawing, transferring and blending a certain volume of liquid. Both sterile and non-sterile Pasteur pipettes are available.

Features

1. NEST Pasteur pipettes are made of LDPE and manufactured in ISO 8 grade cleaning room.
2. Good chemical stability
3. Transparent with clear graduations for easy observation
4. Pipette can be twisted for use with micro vessels and special-shaped containers
5. High precision and good repeatability
6. Individually packed pasteur pipettes are sterilized by EO.
7. DNase/RNase, endotoxin and pyrogenic free

Description

Temperature Range :

Shelf Life :

Product Range

Cat.No.	Volume (mL)	Length (mm)	Pack	Sterile	/Pack	/Box	/Case
318012	1	160	individually packed	Yes	1	500	2,000
318031	1	160	Bagged	No	/	500	5,000

Technical Drawing

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018